


E-Health: Lab2 Software Tools

Robert Martí

robert.marti@udg.edu / D.016 (P4 building)




Aim

- Understand the use of software tools for annotation and visualization.
- Install and use Slicer 3D, and ITK-Snap
- With ITK-Snap
 - View 3D volume (Knee MRI) and ground truth (GT)
 - Perform semi-automatic segmentation.
- With Slicer 3D
 - Evaluate the similarity of the two segmentations using Dice similarity.
- Software:
 - ITK Snap. <http://www.itksnap.org>
 - Slicer 3D. <https://www.slicer.org/>

Itk-Snap


1. Install Itk-Snap & download the 2 knee MRIs.
2. Load the MRI image and segmentation (labels) for image-80.
3. Show the 3D rendering of the labels.



4. Segment the image-81 (only the tibia) with two semi-automatic methods using active contours (clustering, classification or edge attraction).
 1. Save the segmentations into separate files (seg1 & seg2) (use *mha* format).

Slicer

5. Install Slicer
6. Download the extension *DiceComputation* (View-ExtensionsManager)



7. Load (*Add-Data*) the previous segmentations done with Snap. S1 & S2.
8. Using *Volumes*, convert them to labels (*convert to scalar value*)
9. Use the *DiceComputation* to compute the Dice similarity.

Results		
	1	2
1		0.975
2	0.975	

What to submit

- Write down a short summary (5 pages aprox) of your work, including:
 - Problems encountered and snapshots of the correct execution of points 3 (rendering), 4 (segmentation) & 9 (Dice coefficient).
 - Files of the manual segmentations performed